



Vet X

Optimize your Vet Life

Better diagnostic, better treatment

PX 20 HF Plus BT

C'EST L'APPAREIL LE PLUS COMPACT

SON TRÈS FAIBLE POIDS EST UNE CARACTÉRISTIQUE ESSENTIELLE POUR RÉDUIRE LA CHARGE DE TRAVAIL PHYSIQUE LORS DES EXAMENS AVANT L'ACHAT.

**APPAREIL « TOUT TERRAIN »
SÉCURITÉ
POUR LES EXAMENS D'ACHAT ÉQUINS
DOS ET GRASSETS AP
90KV - 20MA
TOTALEMENT CHARGÉ EN 4H
DÉMARRAGE RAPIDE
BATTERIE LITHIUM-POLYMÈRES**



XV192025



**JUSQU'À 560 RADIOS
ROBUSTE ET COMPACT
POUR CR, DR OU ARGENTIQUE
GÉNÉRATEUR PORTABLE
STATIF PAS INDISPENSABLE
MÉMOIRES ET PRESETS
PAS DE CÂBLES
QUALITÉ EXAMION**

RADIOLOGIE PORTABLE



Technical Specifications :

HIGH-FREQUENCY X-RAY GENERATOR :

Tube : Toshiba OX/70 CEI
Output : 90 kV / 20 mA
Weight : 7.2 kg
Adjustable kV : 50 - 90 kV - (1 kV steps)
mAs range : 0.4 - 20 mAs (25 steps)
Current range mA output : 15 - 20 mA
20mA @ 50 ~ 80kV
15mA @ 81 ~ 90kV

Focal spot : 1.2 x 1.2 mm
Total filtration : 2.5 mm Al (including collimator)
Light: LED
Collimator light field : Max : 35 x 35 cm @66 cm SID / 40 x 40 cm@75 cm SID;
Min : 5 x 5 cm @ 100cm SID

Collimator lamp Duration : 30 sec
SID measurement : dual laser pointer & tape measure
Dimensions (L x W x H) mm : 307 x 203 x 174 mm

Battery : Type Li-polymer
Capacity (mAh) 4800
Exposures (approx.) approx. 560
On low exposure volume up to one month per charge
Charging time Less than 4 hours

Included accessories : Memory exposures (15 APR standard)
Aluminum carrying case with shoulder strap

Hybrid : power cable possible
Podoblok/MXR accessory compatible

Warranty : 12 months

