



Better diagnostic, better treatment

# *X-R Static UNIVET*

FOR THE MOST DEMANDING VETERINARIANS

## X-Ray generator

- Table with floating top of 160 cm
- High frequency generator
- Compact and modern design
- Powered by capacitors: a wall outlet is sufficient
- Dual focus for optimal resolution
- Color touch screen control panel



STATIONARY X-RAY



## Technical Specifications :

	X-R Static UNIVET 200HF	X-R Static UNIVET 400HF
Output power	20kW	32kW
Kv	from 40 to 100 kV in steps of 1kV	from 40 to 125 kV in steps of 1kV
ma	up to 300mA	up to 400mA
Tube	0.8/1.8mm	0.6/1.3mm
control screen	15.6 inch color touchscreen	15.6 inch color touchscreen
Table top surface	160x68cm electro-magnetic brakes	160x68cm electro-magnetic brakes
ref.	XV192105	XV192103

