



Vet X
Optimize your Vet Life

Better diagnostic, better treatment

PX 20 HF BT Lite

IT IS A PERFECT CHOICE FOR CONTINUOUS OPERATION

**UNLEASHED ALL-ROUNDER!
SAFETY
EQUINE PRE-PURCHASE
INCLUDING PA STIFLES,
CLEAN QUICK START TO X-RAY
LITHIUM-ION POLYMER**




- Hybrid: the unit can be operated by internal battery or external power supply.
- Great advantage for outdoor radiography
- Innovative high-performance lithium-ion polymer battery approach
- Up to 560 exposures on a single charge
- Clean diagnostic images thanks to high-frequency technology
- Exclusive remote control using a manual switch
- Full charge only within 4 hours

PORTABLE X - RAY



Technical Specifications : PX 20 HF BT Lite

Output	
Tube Vol./Current	90kV / 20 mA
Voltage range	50 -90 kV in 1 kV steps
mAs Range	0,4mAs -20mAs 25 steps
Max. output	1,6 kW at 80 kV
X-ray Tube	
Application	OX / 70 CEI
Focal Spot	1,2mm
HU capacity	10 kHU
Collimator	
	Max. Light field: 35 cm x 35 cm Min. Light field: 5 cm x 5 cm @ 1
Lamp	LED
Laser Pointer	built-in dual laser
Power	
Input	110 / 220 VAC (Free-Voltage)
Frequency/phase	50 / 60 Hz, single phase
Voltage comp.	± 10%
Features	
	7 - segment LED read out reversible, Hand switch
WIFI DR Interface	optional
PROM Memory	15 APR
Weight	7,2 kg
Dimensions	203 x 174 x 307 mm
Battery	
Type	Li-polymer
Capacity (mAh)	4800
Exposures full charge	max. 560 at 80 kV / 2 mAs
Charging time	4 hours



Standard Accessory

Charger, Battery & Hand switch



Aluminium Carrying Case

