

X-DRS VET Expert

Nieuw compleet radiologiesysteem, ideaal voor uw bedrijf

Het radiografiesysteem integreert alle onderdelen die nodig zijn om radiografische beelden te produceren in een compacte eenheid.

- Het radiografieproces wordt door de gebruiker bestuurd via de **all-in-one pc**, zonder dat hij zich hoeft te concentreren op het dier.
- De **bekabeling is volledig weggewerkt** in het systeem, voor maximale bescherming.
- De **EXAMION X-DR XL detector** maakt gebruik van onze **gebruiksvriendelijke X-AQS** radiologiesoftware om briljante röntgenbeelden te genereren voor een duidelijke diagnose.
- **X-AI kunstmatige intelligentie** kan worden geïntegreerd voor extra diagnostische ondersteuning.



Technische kenmerken

Radiosysteem

Dimensions with table top (WxHxD)	160 x 200 x 93 cm
Dimensions with stand (WxHxD)	230 x 200 x 135,3 cm
Weight of the entire system	250 kg
Guaranty	5 years

Software

X-AQS PACS VET	X-AQS PACS VetX-AQS Module Archiving on CD/DVDX-AQS Module Data Compressionincludes for clientX-AQS Module Patient CD/DVDX-AQS Module Import Bitmap ImagesX-AQS Module Import DICOM ImagesX-AQS Module Import AQS ImagesX-AQS Module Image ExportX-AQS Module Mail Export
X-AQS Client	1 workstation (expandable up to 4 workstations)X-AQS ViewerX-AQS DiagnosisLicense per Client, floating licensefor X-AQS Software

Detector

Detector choice	X-DR XL (43 x 43 cm / 17 x 17 In)
Pixel pitch	139 µm
Pixel matrix	3.072 px x 3.072 px
Weight	4,0 kg
Safety class	IPX1
Surface load	150 kg
Point load	100 kg bei Ø 40 mm
Warranty	5 years; Extended warranty on request

Generator

Power	30 kW
Anode	Rotating, max. current 425 mA
Focus	0,6/1,3 mm
Power supply	230 V; 16 A (3 kW)
Ma-range in steps of 25	25 bis 300
KV range in steps of 1KV	40 bis 120
Collimator	R104

